For Research Use Only. Not for use in diagnostic procedures. Product is manufactured in accordance with current Good Manufacturing Practices.

Identification of the reagent

Name	CyLyse™ FX
REF No.	AX032552

Product Description

CyLyse™ FX is a lysing solution with a fixative for red blood cell lysis in the preparation of biological samples from human peripheral blood after staining leukocytes with fluorochrome-conjugated antibodies prior to the flow cytometry analysis.

Principle of the procedure

Leukocyte analysis and detection in peripheral blood requires elimination of interfering cells, mainly erythrocytes. Direct blood sample staining followed by red blood cell lysis is a fast and easy method for whole blood flow cytometry analysis.

Storage and shelf life

Unopened product

Store CyLyse™ FX (10x concentrate) at 2-25 °C. Do not freeze. Do not use after the expiration date stated on the label.

After first opening

CyLyse™ FX retains its performance characteristics after having been placed into use for 36 months.

Components

CyLyse™ FX (10x concentrate) is provided as 50 mL of a proprietary buffered solution containing 20-30 % (v/v) diethylene glycol, < 15 % (v/v) formaldehyde and < 5 % (v/v) methanol.

Evidence of deterioration

Avoid contamination of reagents. In case of component deterioration seen as a visible precipitation or discoloration of the reagent or if data obtained show any performance alteration, please contact the Technical Support of your local Sysmex representative.

Any problem that has occurred in relation to the product shall be reported by the user to the manufacturer.

Precaution and warnings

Important information regarding the safe handling, transport, and disposal of this product is contained in the Safety Data Sheet (available at <http://www.sysmex-partec.com/services>).

Always meet the national and international guidelines and regulatory standards for PPE (personal protective equipment).

Warning symbols

	GHS07		GHS08
--	-------	--	-------

Signal word Danger	
Warnings	
H302	Harmful if swallowed.
H315	Causes skin irritation.
H317	May cause an allergic skin reaction.
H319	Causes serious eye irritation.
H335	May cause respiratory irritation.
H341	Suspected of causing genetic defects.
H350	May cause cancer.
H371	May cause damage to organs.
H373	May cause damage to organs through prolonged or repeated exposure.
Precautions	
P201	Obtain special instructions before use.
P260	Do not breathe dust/fume/gas/mist/vapors/spray.
P280	Wear protective gloves/protective clothing/eye protection/face protection.
P308+P313	IF exposed or concerned: Get medical advice/attention.

Additional required equipment	
Instrument	Flow cytometer equipped with appropriate computer hardware and software. The flow cytometer must be equipped to detect forward scatter (FSC) and side scatter (SSC).
Required laboratory equipment	Vortex mixer, Centrifuge (including rotor system designed for up to 500 g), Material necessary for the collection of whole blood, Disposable test tubes (12 x 75 mm) for staining of samples, Pipettes with disposable tips for 10, 100 and 1000 μL, Adequate personal protective equipment
Required reagents	Fluorochrome-conjugated antibody reagents (e.g. Sysmex CyFlow™ antibody reagents), Phosphate buffered saline (PBS; pH 7.4), Deionized water

Other materials can be required. Refer to the appropriate antibody reagent Instructions for Use (IFU) for more information.

Reagent preparation	
Dilute CyLyse™ FX (10x concentrate) with deionized water at room temperature (1 volume of concentrated solution with 9 volumes of deionized water). Diluted solution is stable for one month when stored at 18-25 °C.	
Disposal	
All disposables, which have been in contact with biohazardous material, must be decontaminated and disposed of according to local legislations and laws. Clean and disinfect contaminated surfaces immediately, use appropriate procedures of decontamination. Always dispose blood samples, assays and accessory fluids after expiration of the maximal storage time.	
Primary sample collection, handling and storage	
WARNING Consider all biological specimens and materials which come in contact with them as biohazardous. Specimens should be handled as potentially infectious and disposed in accordance with federal, state and local regulations.	
Collect whole blood in a sterile tube with K3 or K2 EDTA as anticoagulant. Follow the antibody reagent IFU for sample handling and storage.	
Examination procedure	
Lyse/No-wash procedure	
1. Stain whole blood samples following instructions in the antibody reagent IFU.	
2. Add 500 μL of 10x diluted CyLyse™ FX per 50 μL of specimen and vortex gently.	
3. Incubate for 10-15 minutes at 18-28 °C in the dark.	
4. Analyze sample immediately using flow cytometer or follow instructions in the antibody reagent IFU.	
Lyse/Wash procedure	
1. Stain whole blood samples following instructions in the antibody reagent IFU.	
2. Add 1.0 mL of 10x diluted CyLyse™ FX per 50 μL of specimen and vortex gently.	
3. Incubate for 10-15 minutes at 18-28 °C in the dark.	
4. Centrifuge tubes for 5 minutes at 300 g and remove the supernatant by decanting.	
5. For subsequent analysis on a flow cytometer, resuspend the cell pellet in a sufficient volume of PBS ap-	

propriate for your flow cytometer. For later analysis, follow instructions in the antibody reagent IFU.

Limitations

The test is intended for professional and appropriately trained users in laboratories performing flow cytometry analysis. In case of hyperleukocytosis, it is recommended to dilute blood samples with PBS to a concentration of 5 x 10⁶ leukocytes/mL [1, 2, 3]. In certain disease states, such as haemoglobinopathies, lysis of red cells may be slow, incomplete or even impossible. In this case, it is recommended to isolate mononucleated cells using a density gradient (Ficoll, for example) prior to staining [4, 5, 6, 7, 8, 9]. Samples with nucleated red blood cells may show incomplete lysis of red blood cells. This may also occur when assaying blood samples from patients with certain hematologic disorders in which red cells are difficult to lyse, as in myelofibrosis, sickle-cell anemia or thalassemia [4, 5, 7, 10, 11]. The result for the Lyse/No-wash procedure may vary depending on the analysis platform. It is recommended to validate the Lyse/No-wash procedure for your application and flow cytometer. The flow cytometer may produce false results, if the device has not been aligned and maintained appropriately. Data may be incorrectly interpreted, if fluorescent signals were compensated wrongly or if gates were positioned inaccurately. Accurate and reproducible results will be obtained as long as the procedures used are in accordance with the Technical Data Sheet and compatible with good laboratory practices.


Properly interpreted, if fluorescent signals were compensated wrongly or if gates were positioned inaccurately. Accurate and reproducible results will be obtained as long as the procedures used are in accordance with the Technical Data Sheet and compatible with good laboratory practices.

Literature references

- N. Abramson, B. Melton. Leukocytosis: Basics of Clinical Assessment. American Family Physician, 62 (9): 2053-2060 (2000)
- Kurec. Lipemia and hyperleukocytosis can lead to CBC errors. Medical Laboratory Observer (MLO), 48 (3): 44 (2016)
- L. K. Riley, J. Rupert. Evaluation of Patients with Leukocytosis, American Family Physician, 92 (11): 1004-1011 (2015)
- T. Constantino. Nucleated RBCs—Significance in the Peripheral Blood Film, Laboratory Medicine, Volume 31, Issue 4, 223–229 (2000)
- S. Buoro. Evaluation of nucleated red blood cell count by Sysmex XE-2100 in patients with thalassaemia or sickle cell anaemia and in neonates, Blood Transfus., 13(4): 588-94 (2015)
- F. Booth. Resistance to lysis of erythrocytes containing haemoglobin C—detected in a differential white cell counting system. J Clin Pathol., 36(7): 816–818 (1983)
- Posteraro. The Diagnostic Significance of a Prolonged Erythrocytic Glycerol Lysis Time (GLT50), American Journal of Clinical Pathology, Volume 70, Issue 4, Pages 637–641 (1978)
- E. Genuardi. Ficoll - hypaque separation vs whole blood lysis: Comparison of efficiency and impact on minimal residual disease analysis, Int J Lab Hematol., 40(2):201-208 (2018)
- P. Dagur. Collection, Storage, and Preparation of Human Blood Cells, Curr Protoc Cytom., 73: 5.1.1–5.1.16. (2015)
- S. Buoro. Which clinical significance has automatic detection of very low levels of nucleated red blood cells in the peripheral blood? Ann Transl Med., 4(11): 230 (2016)
- P. Danise. Nucleated red blood cells and soluble transferrin receptor in thalassemia syndromes: relationship with global and ineffective erythropoiesis, Clin Chem Lab Med, 47:1539–42 (2009)

Representative data

The following representative data was obtained using human peripheral whole blood stained with Sysmex Partec CyFlow™ antibody reagents (CD3-APC, CD19-PE and CD45-Pacific Blue™) and treated with CyLyse™ FX. The data was collected on a Sysmex Partec flow cytometer equipped with violet (405 nm), blue (488 nm) and red (638 nm) lasers.

Manufacturer

Sysmex Partec GmbH Am Flugplatz 13 02828 Görlitz Tel.: +49 3581 8746 0 Fax: +49 3581 8746 70 E-mail: info@sysmex-partec.com www.sysmex-partec.com	Distributed in the U.S.A. by: Sysmex America, Inc. 577 Aptakisic Road Lincolnshire, IL 60069, U.S.A. Tel.: 1-888-879-7639 or 1-847-367-3503 E-mail: FCM_Applications@sysmex.com www.sysmex.com/us
--	---

Symbols

	Reference Number		Batch Code		Concentrated Reagent
--	------------------	--	------------	--	----------------------

	Caution		Temperature Limitation		Use By
--	---------	--	------------------------	--	--------

	Legal Manufacturer
--	--------------------

Revision and date of issue	
Rev.:	001
Rev. Date:	03.06.2019

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