CyFlow™ CLEAN

EN

Identification of the IVD reagent

CyFlow™ CLEAN

Intended use

For in vitro diagnostic use only
CyFlow CLEAN is a strong alkaline detergent intended for use as a cleaning agent for the fluidics of Sysmex Flow Cytometry. It is intended to be used for in vitro diagnostic purposes by healthcare professionals and properly trained personnel in a laboratory environment.

Components

Sodium Hypochlorite (active chlorine concentration 5.0 %)

Warnings and precautions

- Follow the warnings and precautions, as described on the product container, package box, package insert or device's Instructions for Use, to handle the product correctly.
 - Do not use the product if it shows signs of contamination or instability, such as turbidity or discoloration.
 - Be sure to wear adequate personal protective equipment, such as protective gloves, a protective mask, protective eyewear, and a lab coat.
 - Avoid direct contact with skin, eyes and mucous membranes. Do not ingest.
 - In case of skin contact, rinse immediately with plenty of water. In case of contact with eyes or mucous membranes, rinse immediately with plenty of water, and seek medical attention. In case of ingestion, seek medical attention immediately.
 - Please see the Safety Data Sheet for the country-specific classification.
- Danger
- | | |
|------|--|
| H314 | Causes severe skin burns and eye damage. |
| P264 | Wash skin thoroughly after handling. |
| P280 | Wear protective gloves/protective clothing/eye protection/face protection. |
- 
- | | |
|----------------|--|
| P303+P361+P353 | IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. |
| P305+P351+P338 | IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. |
| P310 | Immediately call a POISON CENTER/doctor. |
| P501 | Dispose of contents/container to a facility in accordance with local and national regulations. |
| EUH031 | Contact with acids liberates toxic gas. |

Examination procedure

Refer to the device's Instructions for Use.

Storage and shelf life of unopened product

Store the product at 1-30 °C. Keep away from humidity, direct sunlight and heat. When the product is properly stored in its sealed container, it is stable until the expiration date printed on the label.

Storage and shelf life after first opening

Once opened the product is stable for 60 days at 1-30 °C. After opening the product container, store in the box and avoid direct contact with sunlight.

Limitations of the examination procedure

- Always use the product within the specified analysis environment for the device. Accurate analysis results may not be obtained in an environment that does not meet the specified conditions. Refer to the device's Instructions for Use for further information.
- Results may be inaccurate due to interfering substances such as dust, dirt or bacteria growth in the sample or product.
- The reliability of the analysis values cannot be guaranteed if the product is used outside the prescribed intended use.
- The product must not be used after its expiration date.
- Do not refill and reuse product containers.
- Avoid contamination with dust or bacteria after the container is opened. If the product displays any signs of contamination or instability, as indicated by cloudiness or color change, it should be replaced.
- Do not use a product that is suspected to have been frozen.
- The use of this product is validated on specific devices to optimize product performance and meet product specifications. Please refer to the Instructions for Use of your device to confirm that the use of this product is authorized by Sysmex. Sysmex cannot take the responsibility for patient results obtained from the use of Sysmex products on unauthorized devices. It is the responsibility of the user to validate modifications to these instructions or use of the product on devices other than those specified by Sysmex.

Disposal procedures

Disposal procedures for emptied containers, remaining liquids and waste effluents from the device should meet the requirements of applicable local regulations.

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