

Technical Data Sheet

Product

**7-AAD Viability Staining Solution**

Cat. Number/Size

**EXB0026**

**400 tests**

**For Research Use Only.**

**Not for use in diagnostic or therapeutic procedures.**

Application

FC (QC tested)

Description

The 7-AAD (7-aminoactinomycin D) Staining Solution is intended to evaluate cell viability by flow cytometry.

Excitation laser

blue (488 nm)

Storage and handling

Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Images and References

[www.exbio.cz](http://www.exbio.cz)

Unless indicated otherwise, all products are For Research Use Only and not for diagnostic use. In vivo diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s. © 2019 EXBIO Praha, a.s. All rights reserved.